

specimen	tissue	stage	EF1a	optix	dpp	spalt
45	forewing	4	7.08119933	2.61113065	4.78404612	3.52658928
63	forewing	4	7.53959732	2.46389491	5.19402079	3.49436274
66	forewing	4	7.18406609	2.80922007	4.68171628	3.53816727
67	forewing	4	6.38238009	2.09427232	3.97378960	2.78311995
212	forewing	5	6.72314452	2.15559548	4.24710324	2.7989256
216	forewing	5	3.88652651	1.88941700	2.31636374	2.60946775
217	forewing	5	6.98092905	2.31474567	4.54953203	2.89214649
218	forewing	5	6.31832029	1.56820110	3.67695032	1.88315859
45	hindwing	4	7.15207026	1.77754100	4.79332730	3.43818026
63	hindwing	4	7.40649324	2.98186035	5.0810706	3.58869336
66	hindwing	4	7.12332708	2.39207592	4.75589063	3.66513619
67	hindwing	4	6.96529265	2.06901009	4.57764117	2.82535180
212	hindwing	5	6.46965683	2.67525680	4.03297436	2.83067528
216	hindwing	5	6.46521760	1.90634038	3.77635606	2.87559141
217	hindwing	5	6.38266996	2.11657409	4.04261255	2.80364544
218	hindwing	5	5.08392914	1.41899515	2.42953923	1.71578292

Well	Fluor	Target	Sample	tissue	stage	Cq	mean Cq	range	SQ	SQ	Cq	slope
A01	Cy5	EF1a	E6			18.57	17.36	1.83	6.38312102	6.30102999	18.57	-3.29115826
A02	Cy5	EF1a	E6			16.77				6.30102999	16.77	
A03	Cy5	EF1a	E6			16.75				6.30102999	16.75	intercept
A04	Cy5	EF1a	E5			24.44	21.55	4.71	5.11007240	5.30102999	24.44	38.3691786
A05	Cy5	EF1a	E5			20.49				5.30102999	20.49	
A06	Cy5	EF1a	E5			19.73				5.30102999	19.73	Efficiency
A07	Cy5	EF1a	E4			24.91	23.77	1.99	4.43667209	4.30102999	24.91	1.00650145
A08	Cy5	EF1a	E4			23.47				4.30102999	23.47	101%
A09	Cy5	EF1a	E4			22.92				4.30102999	22.92	
B01	Cy5	EF1a	66RF	forewing	4	15.03	14.73	OK	7.18406609	3.30102999	27.64	R^2
B02	Cy5	EF1a				14.56				3.30102999	27.67	0.89447439
B03	Cy5	EF1a				14.59				3.30102999	27.47	
B04	Cy5	EF1a	63RF	forewing	4	13.71	13.56	OK	7.53959732			
B05	Cy5	EF1a				13.42						
B06	Cy5	EF1a				13.54						
B07	Cy5	EF1a	216RF	forewing	5	26.97	25.58	2.48	3.88652651			
B08	Cy5	EF1a				25.28						
B09	Cy5	EF1a				24.49						
C01	Cy5	EF1a	66RH	hindwing	4	15.59	14.93	OK	7.12332708			
C02	Cy5	EF1a				14.67						
C03	Cy5	EF1a				14.51						
C04	Cy5	EF1a	63RF	hindwing	4	15.13	13.99	1.92	7.40649324			
C05	Cy5	EF1a				13.64						
C06	Cy5	EF1a				13.21						
C07	Cy5	EF1a	216RF	hindwing	5	18.41	17.09	2.38	6.46521760			
C08	Cy5	EF1a				16.83						
C09	Cy5	EF1a				16.03						
D01	Cy5	EF1a	67RF	forewing	4	18.06	17.36	OK	6.38238009			
D02	Cy5	EF1a				17.16						
D03	Cy5	EF1a				16.87						
D04	Cy5	EF1a	217RF	forewing	5	17.71	15.39	4.08	6.98092905			

Well	Fluor	Target	Sample	tissue	stage	Cq	mean Cq	range	SQ	SQ	Cq	slope			
D05	Cy5	EF1a	212RF	forewing	5	14.84	16.24	2.26	6.72314452						
D06	Cy5	EF1a				13.63									
D07	Cy5	EF1a				17.40									
D08	Cy5	EF1a				16.19									
D09	Cy5	EF1a	67RF	hindwing	4	15.13	15.45	OK	6.96529265						
E01	Cy5	EF1a				15.50									
E02	Cy5	EF1a				15.39									
E03	Cy5	EF1a				15.44									
E04	Cy5	EF1a	217RF	hindwing	5	20.54	17.36	5.10	6.38266996						
E05	Cy5	EF1a				16.10									
E06	Cy5	EF1a				15.44									
E07	Cy5	EF1a				19.01									
E08	Cy5	EF1a	212RH	hindwing	5	17.01	17.08	3.80	6.46965683						
E09	Cy5	EF1a				15.21									
F01	Cy5	EF1a				15.65									
F02	Cy5	EF1a				14.79									
F03	Cy5	EF1a	45RF	forewing	4	14.76	15.06	OK	7.08119933						
F04	Cy5	EF1a				20.67									
F05	Cy5	EF1a				16.13									
F06	Cy5	EF1a				15.92									
F07	Cy5	EF1a	NRT			35.58	33.21	3.67	1.56737823						
F08	Cy5	EF1a	NRT			32.14									
F09	Cy5	EF1a	NRT			31.91									
G01	Cy5	EF1a	45RH			hindwing							4	15.30	14.83
G02	Cy5	EF1a	14.56												
G03	Cy5	EF1a	14.63												
G04	Cy5	EF1a	218RH	hindwing	5	24.28									
G05	Cy5	EF1a				20.63									
G06	Cy5	EF1a				20.00									
G07	Cy5	EF1a				NTC							#DIV/0!	OK	#DIV/0!
G08	Cy5	EF1a				NTC									

Well	Fluor	Target	Sample	tissue	stage	Cq	mean Cq	range	SQ	SQ	Cq	slope
G09	Cy5	EF1a	NTC									
H01	Cy5	EF1a	E3			27.64	27.59	OK	3.27425446			
H02	Cy5	EF1a	E3			27.67						
H03	Cy5	EF1a	E3			27.47				SQ	Cq	slope
A01	FAM	optix	E6			19.13	17.89	1.91	6.37126139	6.30102999	19.13	-3.34633938
A02	FAM	optix	E6			17.33				6.30102999	17.33	
A03	FAM	optix	E6			17.22				6.30102999	17.22	intercept
A04	FAM	optix	E5			25.04	22.08	4.80	5.12062317	5.30102999	25.04	39.21300861
A05	FAM	optix	E5			20.95				5.30102999	20.95	
A06	FAM	optix	E5			20.24				5.30102999	20.24	Efficiency
A07	FAM	optix	E4			25.45	24.32	2.03	4.45114945	4.30102999	25.45	0.99495630
A08	FAM	optix	E4			24.08				4.30102999	24.08	99%
A09	FAM	optix	E4			23.42				4.30102999	23.42	
B01	FAM	optix	66RF	forewing	4	30.20	29.81	OK	2.80922007	3.30102999	28.21	R^2
B02	FAM	optix				29.58				3.30102999	28.30	0.89381935
B03	FAM	optix				29.65				3.30102999	28.39	
B04	FAM	optix	63RF	forewing	4	31.18	30.97	OK	2.46389491			
B05	FAM	optix				30.81						
B06	FAM	optix				30.92						
B07	FAM	optix	216RF	forewing	5	34.04	32.89	1.97	1.88941700			
B08	FAM	optix				32.56						
B09	FAM	optix				32.07						
C01	FAM	optix	66RH	hindwing	4	31.71	31.21	OK	2.39207592			
C02	FAM	optix				31.06						
C03	FAM	optix				30.85						
C04	FAM	optix	63RF	hindwing	4	30.18	29.23	1.69	2.98186035			
C05	FAM	optix				29.03						
C06	FAM	optix				28.49						
C07	FAM	optix	216RF	hindwing	5	34.14	32.83	2.39	1.90634038			
C08	FAM	optix				32.61						
C09	FAM	optix				31.75						

Well	Fluor	Target	Sample	tissue	stage	Cq	mean Cq	range	SQ	SQ	Cq	slope
D01	FAM	optix	67RF	forewing	4	32.78	32.20	OK	2.09427232i			
D02	FAM	optix				31.77						
D03	FAM	optix				32.07						
D04	FAM	optix	217RF	forewing	5	33.57	31.47	3.71	2.31474567i			
D05	FAM	optix				30.97						
D06	FAM	optix				29.86						
D07	FAM	optix	212RF	forewing	5	33.09	32.00	1.97	2.15559548i			
D08	FAM	optix				31.78						
D09	FAM	optix				31.12						
E01	FAM	optix	67RF	hindwing	4	33.28	32.29	1.63	2.06901009			
E02	FAM	optix				31.93						
E03	FAM	optix				31.65						
E04	FAM	optix	217RF	hindwing	5	34.99	32.13	4.53	2.11657409i			
E05	FAM	optix				30.94						
E06	FAM	optix				30.46						
E07	FAM	optix	212RH	hindwing	5	32.19	30.26	3.71	2.67525680i			
E08	FAM	optix				30.11						
E09	FAM	optix				28.48						
F01	FAM	optix	45RF	forewing	4	31.01	30.48	OK	2.61113065i			
F02	FAM	optix				30.01						
F03	FAM	optix				30.41						
F04	FAM	optix	218RF	forewing	5		33.97	OK	1.56820110i			
F05	FAM	optix				34.33						
F06	FAM	optix				33.60						
F07	FAM	optix	NRT			36.27	33.83	4.20	1.60760334i			
F08	FAM	optix	NRT			33.17						
F09	FAM	optix	NRT			32.07						
G01	FAM	optix	45RH	hindwing	4	33.67	33.26	OK	1.77754100i			
G02	FAM	optix				33.27						
G03	FAM	optix				32.86						
G04	FAM	optix	218RH	hindwing	5	38.27	34.46	5.90	1.41899515			

Well	Fluor	Target	Sample	tissue	stage	Cq	mean Cq	range	SQ	SQ	Cq	slope
G05	FAM	optix				32.76						
G06	FAM	optix				32.37						
G07	FAM	optix	NTC				#DIV/0!	OK	#DIV/0!			
G08	FAM	optix	NTC									
G09	FAM	optix	NTC									
H01	FAM	optix	E3			28.21	28.30	OK	3.26108596			
H02	FAM	optix	E3			28.30						
H03	FAM	optix	E3			28.39				SQ	Cq	slope
A01	HEX	spalt	E6			19.41	18.20	1.85	6.38313052	6.30102999	19.41	-3.35032625
A02	HEX	spalt	E6			17.62				6.30102999	17.62	
A03	HEX	spalt	E6			17.57				6.30102999	17.57	intercept
A04	HEX	spalt	E5			25.40	22.48	4.77	5.10578844	5.30102999	25.40	39.5870685
A05	HEX	spalt	E5			21.41				5.30102999	21.41	
A06	HEX	spalt	E5			20.63				5.30102999	20.63	Efficiency
A07	HEX	spalt	E4			25.80	24.69	1.96	4.44521152	4.30102999	25.80	0.99414193
A08	HEX	spalt	E4			24.44				4.30102999	24.44	99%
A09	HEX	spalt	E4			23.84				4.30102999	23.84	
B01	HEX	spalt	66RF	forewing	4	28.01	27.73	OK	3.53816727	3.30102999	28.61	R^2
B02	HEX	spalt				27.61				3.30102999	28.67	0.89577073
B03	HEX	spalt				27.58				3.30102999	28.61	
B04	HEX	spalt	63RF	forewing	4	28.10	27.88	OK	3.49436274			
B05	HEX	spalt				27.77						
B06	HEX	spalt				27.77						
B07	HEX	spalt	216RF	forewing	5	32.17	30.84	2.50	2.60946775			
B08	HEX	spalt				30.69						
B09	HEX	spalt				29.67						
C01	HEX	spalt	66RH	hindwing	4	27.87	27.31	OK	3.66513619			
C02	HEX	spalt				27.11						
C03	HEX	spalt				26.94						
C04	HEX	spalt	63RF	hindwing	4	28.59	27.56	1.80	3.58869336			
C05	HEX	spalt				27.30						

Well	Fluor	Target	Sample	tissue	stage	Cq	mean Cq	range	SQ	SQ	Cq	slope
C06	HEX	spalt				26.80						
C07	HEX	spalt	216RF	hindwing	5	31.28	29.95	2.43	2.87559141			
C08	HEX	spalt				29.73						
C09	HEX	spalt				28.85						
D01	HEX	spalt	67RF	forewing	4	30.89	30.26	OK	2.78311995			
D02	HEX	spalt				30.10						
D03	HEX	spalt				29.79						
D04	HEX	spalt	217RF	forewing	5	32.18	29.90	4.07	2.89214649			
D05	HEX	spalt				29.41						
D06	HEX	spalt				28.11						
D07	HEX	spalt	212RF	forewing	5	31.28	30.21	2.11	2.7989256			
D08	HEX	spalt				30.19						
D09	HEX	spalt				29.17						
E01	HEX	spalt	67RF	hindwing	4	30.25	30.12	OK	2.82535180			
E02	HEX	spalt				29.99						
E03	HEX	spalt				30.13						
E04	HEX	spalt	217RF	hindwing	5	33.45	30.19	5.25	2.80364544			
E05	HEX	spalt				28.94						
E06	HEX	spalt				28.20						
E07	HEX	spalt	212RH	hindwing	5	31.98	30.10	3.75	2.83067528			
E08	HEX	spalt				30.10						
E09	HEX	spalt				28.23						
F01	HEX	spalt	45RF	forewing	4	28.36	27.77	OK	3.52658928			
F02	HEX	spalt				27.46						
F03	HEX	spalt				27.50						
F04	HEX	spalt	218RF	forewing	5	37.26	33.28	6.06	1.88315859			
F05	HEX	spalt				31.37						
F06	HEX	spalt				31.20						
F07	HEX	spalt	NRT			36.07	33.85	3.83	1.71122271			
F08	HEX	spalt	NRT			33.25						
F09	HEX	spalt	NRT			32.24						

Well	Fluor	Target	Sample	tissue	stage	Cq	mean Cq	range	SQ	SQ	Cq	slope
G01	HEX	spalt	45RH	hindwing	4	28.51	28.07	OK	3.43818026			
G02	HEX	spalt				27.82						
G03	HEX	spalt				27.88						
G04	HEX	spalt	218RH	hindwing	5	37.13	33.84	5.30	1.71578292			
G05	HEX	spalt				32.55						
G06	HEX	spalt				31.83						
G07	HEX	spalt	NTC				#DIV/0!	OK	#DIV/0!			
G08	HEX	spalt	NTC									
G09	HEX	spalt	NTC									
H01	HEX	spalt	E3			28.61	28.63	OK	3.26998949			
H02	HEX	spalt	E3			28.67						
H03	HEX	spalt	E3			28.61				SQ	Cq	slope
A01	Texas Red	dpp	E6			19.81	18.59	1.84	6.37964774	6.30102999	19.81	-3.32015375
A02	Texas Red	dpp	E6			17.99				6.30102999	17.99	
A03	Texas Red	dpp	E6			17.97				6.30102999	17.97	intercept
A04	Texas Red	dpp	E5			25.72	22.80	4.78	5.11271310	5.30102999	25.72	39.7733852
A05	Texas Red	dpp	E5			21.73				5.30102999	21.73	
A06	Texas Red	dpp	E5			20.94				5.30102999	20.94	Efficiency
A07	Texas Red	dpp	E4			26.08	25.03	1.88	4.44181054	4.30102999	26.08	1.00037049
A08	Texas Red	dpp	E4			24.80				4.30102999	24.80	100%
A09	Texas Red	dpp	E4			24.20				4.30102999	24.20	
B01	Texas Red	dpp	66RF	forewing	4	24.44	24.23	OK	4.68171628	3.30102999	28.87	R^2
B02	Texas Red	dpp				24.10				3.30102999	29.03	0.89560127
B03	Texas Red	dpp				24.14				3.30102999	28.85	
B04	Texas Red	dpp	63RF	forewing	4	22.70	22.53	OK	5.19402079			
B05	Texas Red	dpp				22.41						
B06	Texas Red	dpp				22.47						
B07	Texas Red	dpp	216RF	forewing	5	33.05	32.08	1.94	2.31636374			
B08	Texas Red	dpp				32.09						
B09	Texas Red	dpp				31.11						
C01	Texas Red	dpp	66RH	hindwing	4	24.58	23.98	OK	4.75589063			

Well	Fluor	Target	Sample	tissue	stage	Cq	mean Cq	range	SQ	SQ	Cq	slope
C02	Texas Red	dpp				23.75						
C03	Texas Red	dpp				23.62						
C04	Texas Red	dpp	63RF	hindwing	4	23.97	22.90	1.83	5.0810706			
C05	Texas Red	dpp				22.61						
C06	Texas Red	dpp				22.13						
C07	Texas Red	dpp	216RF	hindwing	5	28.50	27.24	2.35	3.77635606			
C08	Texas Red	dpp				27.06						
C09	Texas Red	dpp				26.15						
D01	Texas Red	dpp	67RF	forewing	4	27.27	26.58	OK	3.97378960			
D02	Texas Red	dpp				26.32						
D03	Texas Red	dpp				26.14						
D04	Texas Red	dpp	217RF	forewing	5	27.06	24.67	4.23	4.54953203			
D05	Texas Red	dpp				24.12						
D06	Texas Red	dpp				22.83						
D07	Texas Red	dpp	212RF	forewing	5	26.82	25.67	2.27	4.24710324			
D08	Texas Red	dpp				25.65						
D09	Texas Red	dpp				24.55						
E01	Texas Red	dpp	67RF	hindwing	4	24.67	24.57	OK	4.57764117			
E02	Texas Red	dpp				24.49						
E03	Texas Red	dpp				24.56						
E04	Texas Red	dpp	217RF	hindwing	5	29.61	26.35	5.21	4.04261255			
E05	Texas Red	dpp				25.05						
E06	Texas Red	dpp				24.40						
E07	Texas Red	dpp	212RH	hindwing	5	28.30	26.38	3.83	4.03297436			
E08	Texas Red	dpp				26.39						
E09	Texas Red	dpp				24.46						
F01	Texas Red	dpp	45RF	forewing	4	24.54	23.89	OK	4.78404612			
F02	Texas Red	dpp				23.61						
F03	Texas Red	dpp				23.52						
F04	Texas Red	dpp	218RF	forewing	5	30.61	27.57	4.73	3.67695032			
F05	Texas Red	dpp				26.20						

Well	Fluor	Target	Sample	tissue	stage	Cq	mean Cq	range	SQ	SQ	Cq	slope
F06	Texas Red	dpp				25.88						
F07	Texas Red	dpp	NRT			35.90	33.83	3.35	1.78902780			
F08	Texas Red	dpp	NRT			33.05						
F09	Texas Red	dpp	NRT			32.55						
G01	Texas Red	dpp	45RH	hindwing	4	24.35	23.86	OK	4.79332730			
G02	Texas Red	dpp				23.56						
G03	Texas Red	dpp				23.66						
G04	Texas Red	dpp	218RH	hindwing	5	34.29	31.71	4.33	2.42953923			
G05	Texas Red	dpp				30.87						
G06	Texas Red	dpp				29.96						
G07	Texas Red	dpp	NTC				38.53	OK	0.37390432			
G08	Texas Red	dpp	NTC			38.29						
G09	Texas Red	dpp	NTC			38.78						
H01	Texas Red	dpp	E3			28.87	28.92	OK	3.26994859			
H02	Texas Red	dpp	E3			29.03						
H03	Texas Red	dpp	E3			28.85						

:

v

:

✓

:

✓

;

✓

File Name CHF_2025-04-17 15-46-43_EcoDevoCFX96.pcrd

Created By I admin

Notes

ID

Run Started 04/17/2025 19:47:23 UTC

Run Ended 04/17/2025 20:06:29 UTC

Sample Vol 20

Lid Temp 105

Protocol File iQ_Probes_Multiplex.prcf

Plate Setup QuickPlate_96 wells_All Channels.pltd

Base Serial I 795BR20257

Optical Hear 795BR20257

CFX Maestr 5.3.022.1030.

Well group All Wells

Amplification 3

Melt step

Well	Fluor	Target	Content	Sample	Cq
A01	Cy5		Unkn		18.57
A02	Cy5		Unkn		16.77
A03	Cy5		Unkn		16.75
A04	Cy5		Unkn		24.44
A05	Cy5		Unkn		20.49
A06	Cy5		Unkn		19.73
A07	Cy5		Unkn		24.91
A08	Cy5		Unkn		23.47
A09	Cy5		Unkn		22.92
B01	Cy5		Unkn		15.03
B02	Cy5		Unkn		14.56
B03	Cy5		Unkn		14.59

B04	Cy5	Unkn	13.71
B05	Cy5	Unkn	13.42
B06	Cy5	Unkn	13.54
B07	Cy5	Unkn	26.97
B08	Cy5	Unkn	25.28
B09	Cy5	Unkn	24.49
C01	Cy5	Unkn	15.59
C02	Cy5	Unkn	14.67
C03	Cy5	Unkn	14.51
C04	Cy5	Unkn	15.13
C05	Cy5	Unkn	13.64
C06	Cy5	Unkn	13.21
C07	Cy5	Unkn	18.41
C08	Cy5	Unkn	16.83
C09	Cy5	Unkn	16.03
D01	Cy5	Unkn	18.06
D02	Cy5	Unkn	17.16
D03	Cy5	Unkn	16.87
D04	Cy5	Unkn	17.71
D05	Cy5	Unkn	14.84
D06	Cy5	Unkn	13.63
D07	Cy5	Unkn	17.40
D08	Cy5	Unkn	16.19
D09	Cy5	Unkn	15.13
E01	Cy5	Unkn	15.50
E02	Cy5	Unkn	15.39
E03	Cy5	Unkn	15.44
E04	Cy5	Unkn	20.54
E05	Cy5	Unkn	16.10
E06	Cy5	Unkn	15.44
E07	Cy5	Unkn	19.01

E08	Cy5	Unkn	17.01
E09	Cy5	Unkn	15.21
F01	Cy5	Unkn	15.65
F02	Cy5	Unkn	14.79
F03	Cy5	Unkn	14.76
F04	Cy5	Unkn	20.67
F05	Cy5	Unkn	16.13
F06	Cy5	Unkn	15.92
F07	Cy5	Unkn	35.58
F08	Cy5	Unkn	32.14
F09	Cy5	Unkn	31.91
G01	Cy5	Unkn	15.30
G02	Cy5	Unkn	14.56
G03	Cy5	Unkn	14.63
G04	Cy5	Unkn	24.28
G05	Cy5	Unkn	20.63
G06	Cy5	Unkn	20.00
G07	Cy5	Unkn	
G08	Cy5	Unkn	
G09	Cy5	Unkn	
H01	Cy5	Unkn	27.64
H02	Cy5	Unkn	27.67
H03	Cy5	Unkn	27.47
A01	Cy5-5	Unkn	
A02	Cy5-5	Unkn	
A03	Cy5-5	Unkn	
A04	Cy5-5	Unkn	
A05	Cy5-5	Unkn	
A06	Cy5-5	Unkn	
A07	Cy5-5	Unkn	
A08	Cy5-5	Unkn	

A09	Cy5-5	Unkn
B01	Cy5-5	Unkn
B02	Cy5-5	Unkn
B03	Cy5-5	Unkn
B04	Cy5-5	Unkn
B05	Cy5-5	Unkn
B06	Cy5-5	Unkn
B07	Cy5-5	Unkn
B08	Cy5-5	Unkn
B09	Cy5-5	Unkn
C01	Cy5-5	Unkn
C02	Cy5-5	Unkn
C03	Cy5-5	Unkn
C04	Cy5-5	Unkn
C05	Cy5-5	Unkn
C06	Cy5-5	Unkn
C07	Cy5-5	Unkn
C08	Cy5-5	Unkn
C09	Cy5-5	Unkn
D01	Cy5-5	Unkn
D02	Cy5-5	Unkn
D03	Cy5-5	Unkn
D04	Cy5-5	Unkn
D05	Cy5-5	Unkn
D06	Cy5-5	Unkn
D07	Cy5-5	Unkn
D08	Cy5-5	Unkn
D09	Cy5-5	Unkn
E01	Cy5-5	Unkn
E02	Cy5-5	Unkn
E03	Cy5-5	Unkn

E04	Cy5-5	Unkn	
E05	Cy5-5	Unkn	
E06	Cy5-5	Unkn	
E07	Cy5-5	Unkn	
E08	Cy5-5	Unkn	
E09	Cy5-5	Unkn	
F01	Cy5-5	Unkn	
F02	Cy5-5	Unkn	
F03	Cy5-5	Unkn	
F04	Cy5-5	Unkn	
F05	Cy5-5	Unkn	
F06	Cy5-5	Unkn	
F07	Cy5-5	Unkn	
F08	Cy5-5	Unkn	
F09	Cy5-5	Unkn	
G01	Cy5-5	Unkn	
G02	Cy5-5	Unkn	
G03	Cy5-5	Unkn	
G04	Cy5-5	Unkn	
G05	Cy5-5	Unkn	
G06	Cy5-5	Unkn	
G07	Cy5-5	Unkn	
G08	Cy5-5	Unkn	
G09	Cy5-5	Unkn	
H01	Cy5-5	Unkn	
H02	Cy5-5	Unkn	
H03	Cy5-5	Unkn	
A01	FAM	Unkn	19.13
A02	FAM	Unkn	17.33
A03	FAM	Unkn	17.22
A04	FAM	Unkn	25.04

A05	FAM	Unkn	20.95
A06	FAM	Unkn	20.24
A07	FAM	Unkn	25.45
A08	FAM	Unkn	24.08
A09	FAM	Unkn	23.42
B01	FAM	Unkn	30.20
B02	FAM	Unkn	29.58
B03	FAM	Unkn	29.65
B04	FAM	Unkn	31.18
B05	FAM	Unkn	30.81
B06	FAM	Unkn	30.92
B07	FAM	Unkn	34.04
B08	FAM	Unkn	32.56
B09	FAM	Unkn	32.07
C01	FAM	Unkn	31.71
C02	FAM	Unkn	31.06
C03	FAM	Unkn	30.85
C04	FAM	Unkn	30.18
C05	FAM	Unkn	29.03
C06	FAM	Unkn	28.49
C07	FAM	Unkn	34.14
C08	FAM	Unkn	32.61
C09	FAM	Unkn	31.75
D01	FAM	Unkn	32.78
D02	FAM	Unkn	31.77
D03	FAM	Unkn	32.07
D04	FAM	Unkn	33.57
D05	FAM	Unkn	30.97
D06	FAM	Unkn	29.86
D07	FAM	Unkn	33.09
D08	FAM	Unkn	31.78

D09	FAM	Unkn	31.12
E01	FAM	Unkn	33.28
E02	FAM	Unkn	31.93
E03	FAM	Unkn	31.65
E04	FAM	Unkn	34.99
E05	FAM	Unkn	30.94
E06	FAM	Unkn	30.46
E07	FAM	Unkn	32.19
E08	FAM	Unkn	30.11
E09	FAM	Unkn	28.48
F01	FAM	Unkn	31.01
F02	FAM	Unkn	30.01
F03	FAM	Unkn	30.41
F04	FAM	Unkn	
F05	FAM	Unkn	34.33
F06	FAM	Unkn	33.60
F07	FAM	Unkn	36.27
F08	FAM	Unkn	33.17
F09	FAM	Unkn	32.07
G01	FAM	Unkn	33.67
G02	FAM	Unkn	33.27
G03	FAM	Unkn	32.86
G04	FAM	Unkn	38.27
G05	FAM	Unkn	32.76
G06	FAM	Unkn	32.37
G07	FAM	Unkn	
G08	FAM	Unkn	
G09	FAM	Unkn	
H01	FAM	Unkn	28.21
H02	FAM	Unkn	28.30
H03	FAM	Unkn	28.39

A01	HEX	Unkn	19.41
A02	HEX	Unkn	17.62
A03	HEX	Unkn	17.57
A04	HEX	Unkn	25.40
A05	HEX	Unkn	21.41
A06	HEX	Unkn	20.63
A07	HEX	Unkn	25.80
A08	HEX	Unkn	24.44
A09	HEX	Unkn	23.84
B01	HEX	Unkn	28.01
B02	HEX	Unkn	27.61
B03	HEX	Unkn	27.58
B04	HEX	Unkn	28.10
B05	HEX	Unkn	27.77
B06	HEX	Unkn	27.77
B07	HEX	Unkn	32.17
B08	HEX	Unkn	30.69
B09	HEX	Unkn	29.67
C01	HEX	Unkn	27.87
C02	HEX	Unkn	27.11
C03	HEX	Unkn	26.94
C04	HEX	Unkn	28.59
C05	HEX	Unkn	27.30
C06	HEX	Unkn	26.80
C07	HEX	Unkn	31.28
C08	HEX	Unkn	29.73
C09	HEX	Unkn	28.85
D01	HEX	Unkn	30.89
D02	HEX	Unkn	30.10
D03	HEX	Unkn	29.79
D04	HEX	Unkn	32.18

D05	HEX	Unkn	29.41
D06	HEX	Unkn	28.11
D07	HEX	Unkn	31.28
D08	HEX	Unkn	30.19
D09	HEX	Unkn	29.17
E01	HEX	Unkn	30.25
E02	HEX	Unkn	29.99
E03	HEX	Unkn	30.13
E04	HEX	Unkn	33.45
E05	HEX	Unkn	28.94
E06	HEX	Unkn	28.20
E07	HEX	Unkn	31.98
E08	HEX	Unkn	30.10
E09	HEX	Unkn	28.23
F01	HEX	Unkn	28.36
F02	HEX	Unkn	27.46
F03	HEX	Unkn	27.50
F04	HEX	Unkn	37.26
F05	HEX	Unkn	31.37
F06	HEX	Unkn	31.20
F07	HEX	Unkn	36.07
F08	HEX	Unkn	33.25
F09	HEX	Unkn	32.24
G01	HEX	Unkn	28.51
G02	HEX	Unkn	27.82
G03	HEX	Unkn	27.88
G04	HEX	Unkn	37.13
G05	HEX	Unkn	32.55
G06	HEX	Unkn	31.83
G07	HEX	Unkn	
G08	HEX	Unkn	

G09	HEX	Unkn	
H01	HEX	Unkn	28.61
H02	HEX	Unkn	28.67
H03	HEX	Unkn	28.61
A01	Texas Red	Unkn	19.81
A02	Texas Red	Unkn	17.99
A03	Texas Red	Unkn	17.97
A04	Texas Red	Unkn	25.72
A05	Texas Red	Unkn	21.73
A06	Texas Red	Unkn	20.94
A07	Texas Red	Unkn	26.08
A08	Texas Red	Unkn	24.80
A09	Texas Red	Unkn	24.20
B01	Texas Red	Unkn	24.44
B02	Texas Red	Unkn	24.10
B03	Texas Red	Unkn	24.14
B04	Texas Red	Unkn	22.70
B05	Texas Red	Unkn	22.41
B06	Texas Red	Unkn	22.47
B07	Texas Red	Unkn	33.05
B08	Texas Red	Unkn	32.09
B09	Texas Red	Unkn	31.11
C01	Texas Red	Unkn	24.58
C02	Texas Red	Unkn	23.75
C03	Texas Red	Unkn	23.62
C04	Texas Red	Unkn	23.97
C05	Texas Red	Unkn	22.61
C06	Texas Red	Unkn	22.13
C07	Texas Red	Unkn	28.50
C08	Texas Red	Unkn	27.06
C09	Texas Red	Unkn	26.15

D01	Texas Red	Unkn	27.27
D02	Texas Red	Unkn	26.32
D03	Texas Red	Unkn	26.14
D04	Texas Red	Unkn	27.06
D05	Texas Red	Unkn	24.12
D06	Texas Red	Unkn	22.83
D07	Texas Red	Unkn	26.82
D08	Texas Red	Unkn	25.65
D09	Texas Red	Unkn	24.55
E01	Texas Red	Unkn	24.67
E02	Texas Red	Unkn	24.49
E03	Texas Red	Unkn	24.56
E04	Texas Red	Unkn	29.61
E05	Texas Red	Unkn	25.05
E06	Texas Red	Unkn	24.40
E07	Texas Red	Unkn	28.30
E08	Texas Red	Unkn	26.39
E09	Texas Red	Unkn	24.46
F01	Texas Red	Unkn	24.54
F02	Texas Red	Unkn	23.61
F03	Texas Red	Unkn	23.52
F04	Texas Red	Unkn	30.61
F05	Texas Red	Unkn	26.20
F06	Texas Red	Unkn	25.88
F07	Texas Red	Unkn	35.90
F08	Texas Red	Unkn	33.05
F09	Texas Red	Unkn	32.55
G01	Texas Red	Unkn	24.35
G02	Texas Red	Unkn	23.56
G03	Texas Red	Unkn	23.66
G04	Texas Red	Unkn	34.29

G05	Texas Red	Unkn	30.87
G06	Texas Red	Unkn	29.96
G07	Texas Red	Unkn	
G08	Texas Red	Unkn	38.29
G09	Texas Red	Unkn	38.78
H01	Texas Red	Unkn	28.87
H02	Texas Red	Unkn	29.03
H03	Texas Red	Unkn	28.85

Well	Fluor	Target	Sample	tissue	stage	Cq	mean Cq	range	SQ	SQ	Cq	slope
A01	Cy5	EF1a	E6			18.57	17.36	1.83	6.38312102	6.30102999	18.57	-3.29115826
A02	Cy5	EF1a	E6			16.77				6.30102999	16.77	
A03	Cy5	EF1a	E6			16.75				6.30102999	16.75	intercept
A04	Cy5	EF1a	E5			24.44	21.55	4.71	5.11007240	5.30102999	24.44	38.3691786
A05	Cy5	EF1a	E5			20.49				5.30102999	20.49	
A06	Cy5	EF1a	E5			19.73				5.30102999	19.73	Efficiency
A07	Cy5	EF1a	E4			24.91	23.77	1.99	4.43667209	4.30102999	24.91	1.00650145
A08	Cy5	EF1a	E4			23.47				4.30102999	23.47	101%
A09	Cy5	EF1a	E4			22.92				4.30102999	22.92	
B01	Cy5	EF1a	66RF	forewing	4	15.03	14.73	OK	7.18406609	3.30102999	27.64	R^2
B02	Cy5	EF1a				14.56				3.30102999	27.67	0.89447439
B03	Cy5	EF1a				14.59				3.30102999	27.47	
B04	Cy5	EF1a	63RF	forewing	4	13.71	13.56	OK	7.53959732			
B05	Cy5	EF1a				13.42						
B06	Cy5	EF1a				13.54						
C01	Cy5	EF1a	66RH	hindwing	4	15.59	14.93	OK	7.12332708			
C02	Cy5	EF1a				14.67						
C03	Cy5	EF1a				14.51						
C04	Cy5	EF1a	63RH	hindwing	4	15.13	13.99	1.92	7.40649324			
C05	Cy5	EF1a				13.64						
C06	Cy5	EF1a				13.21						
C07	Cy5	EF1a	216RH	hindwing	5	18.41	17.09	2.38	6.46521760			
C08	Cy5	EF1a				16.83						
C09	Cy5	EF1a				16.03						
D01	Cy5	EF1a	67RF	forewing	4	18.06	17.36	OK	6.38238009			
D02	Cy5	EF1a				17.16						
D03	Cy5	EF1a				16.87						
D04	Cy5	EF1a	217RF	forewing	5	17.71	15.39	4.08	6.98092905			
D05	Cy5	EF1a				14.84						
D06	Cy5	EF1a				13.63						
D07	Cy5	EF1a	212RF	forewing	5	17.40	16.24	2.26	6.72314452			

Well	Fluor	Target	Sample	tissue	stage	Cq	mean Cq	range	SQ	SQ	Cq	slope
D08	Cy5	EF1a				16.19						
D09	Cy5	EF1a				15.13						
E01	Cy5	EF1a	67RH	hindwing	4	15.50	15.45	OK	6.96529265			
E02	Cy5	EF1a				15.39						
E03	Cy5	EF1a				15.44						
E04	Cy5	EF1a	217RH	hindwing	5	20.54	17.36	5.10	6.38266996			
E05	Cy5	EF1a				16.10						
E06	Cy5	EF1a				15.44						
E07	Cy5	EF1a	212RH	hindwing	5	19.01	17.08	3.80	6.46965683			
E08	Cy5	EF1a				17.01						
E09	Cy5	EF1a				15.21						
F01	Cy5	EF1a	45RF	forewing	4	15.65	15.06	OK	7.08119933			
F02	Cy5	EF1a				14.79						
F03	Cy5	EF1a				14.76						
F04	Cy5	EF1a	218RF	forewing	5	20.67	17.57	4.75	6.31832029			
F05	Cy5	EF1a				16.13						
F06	Cy5	EF1a				15.92						
F07	Cy5	EF1a	NRT			35.58	33.21	3.67	1.56737823			
F08	Cy5	EF1a	NRT			32.14						
F09	Cy5	EF1a	NRT			31.91						
G01	Cy5	EF1a	45RH	hindwing	4	15.30	14.83	OK	7.15207026			
G02	Cy5	EF1a				14.56						
G03	Cy5	EF1a				14.63						
G04	Cy5	EF1a	218RH	hindwing	5	24.28	21.64	4.27	5.08392914			
G05	Cy5	EF1a				20.63						
G06	Cy5	EF1a				20.00						
G07	Cy5	EF1a	NTC				#DIV/0!	OK	#DIV/0!			
G08	Cy5	EF1a	NTC									
G09	Cy5	EF1a	NTC									
H01	Cy5	EF1a	E3			27.64	27.59	OK	3.27425446			
H02	Cy5	EF1a	E3			27.67						

Well	Fluor	Target	Sample	tissue	stage	Cq	mean Cq	range	SQ	SQ	Cq	slope
H03	Cy5	EF1a	E3			27.47						
A01	FAM	optix	E6			19.13	17.89	1.91	6.37126139			
A02	FAM	optix	E6			17.33						
A03	FAM	optix	E6			17.22				SQ	Cq	slope
A04	FAM	optix	E5			25.04	22.08	4.80	5.12062317	6.30102999	19.13	-3.34633938
A05	FAM	optix	E5			20.95				6.30102999	17.33	
A06	FAM	optix	E5			20.24				6.30102999	17.22	intercept
A07	FAM	optix	E4			25.45	24.32	2.03	4.45114945	5.30102999	25.04	39.21300861
A08	FAM	optix	E4			24.08				5.30102999	20.95	
A09	FAM	optix	E4			23.42				5.30102999	20.24	Efficiency
B01	FAM	optix	66RF	forewing	4	30.20	29.81	OK	2.80922007	4.30102999	25.45	0.99495630
B02	FAM	optix				29.58				4.30102999	24.08	99%
B03	FAM	optix				29.65				4.30102999	23.42	
B04	FAM	optix	63RF	forewing	4	31.18	30.97	OK	2.46389491	3.30102999	28.21	R^2
B05	FAM	optix				30.81				3.30102999	28.30	0.89381935
B06	FAM	optix				30.92				3.30102999	28.39	
B07	FAM	optix	216RF	forewing	5	34.04	32.89	1.97	1.88941700			
B08	FAM	optix				32.56						
B09	FAM	optix				32.07						
C01	FAM	optix	66RH	hindwing	4	31.71	31.21	OK	2.39207592			
C02	FAM	optix				31.06						
C03	FAM	optix				30.85						
C04	FAM	optix	63RF	hindwing	4	30.18	29.23	1.69	2.98186035			
C05	FAM	optix				29.03						
C06	FAM	optix				28.49						
C07	FAM	optix	216RH	hindwing	5	34.14	32.83	2.39	1.90634038			
C08	FAM	optix				32.61						
C09	FAM	optix				31.75						
D01	FAM	optix	67RF	forewing	4	32.78	32.20	OK	2.09427232			
D02	FAM	optix				31.77						
D03	FAM	optix				32.07						

Well	Fluor	Target	Sample	tissue	stage	Cq	mean Cq	range	SQ	SQ	Cq	slope
D04	FAM	optix	217RF	forewing	5	33.57	31.47	3.71	2.314745671			
D05	FAM	optix				30.97						
D06	FAM	optix				29.86						
D07	FAM	optix	212RF	forewing	5	33.09	32.00	1.97	2.155595481			
D08	FAM	optix				31.78						
D09	FAM	optix				31.12						
E01	FAM	optix	67RH	hindwing	4	33.28	32.29	1.63	2.06901009			
E02	FAM	optix				31.93						
E03	FAM	optix				31.65						
E04	FAM	optix	217RH	hindwing	5	34.99	32.13	4.53	2.116574097			
E05	FAM	optix				30.94						
E06	FAM	optix				30.46						
E07	FAM	optix	212RH	hindwing	5	32.19	30.26	3.71	2.675256801			
E08	FAM	optix				30.11						
E09	FAM	optix				28.48						
F01	FAM	optix	45RF	forewing	4	31.01	30.48	OK	2.611130651			
F02	FAM	optix				30.01						
F03	FAM	optix				30.41						
F04	FAM	optix	218RF	forewing	5		33.97	OK	1.568201101			
F05	FAM	optix				34.33						
F06	FAM	optix				33.60						
F07	FAM	optix	NRT			36.27	33.83	4.20	1.607603341			
F08	FAM	optix	NRT			33.17						
F09	FAM	optix	NRT			32.07						
G01	FAM	optix	45RH	hindwing	4	33.67	33.26	OK	1.777541001			
G02	FAM	optix				33.27						
G03	FAM	optix				32.86						
G04	FAM	optix	218RH	hindwing	5	38.27	34.46	5.90	1.41899515			
G05	FAM	optix				32.76						
G06	FAM	optix				32.37						
G07	FAM	optix	NTC				#DIV/0!	OK	#DIV/0!			

Well	Fluor	Target	Sample	tissue	stage	Cq	mean Cq	range	SQ	SQ	Cq	slope
G08	FAM	optix	NTC									
G09	FAM	optix	NTC									
H01	FAM	optix	E3			28.21	28.30	OK	3.26108596			
H02	FAM	optix	E3			28.30						
H03	FAM	optix	E3			28.39						
A01	HEX	spalt	E6			19.41	18.20	1.85	6.38313052			
A02	HEX	spalt	E6			17.62						
A03	HEX	spalt	E6			17.57				SQ	Cq	slope
A04	HEX	spalt	E5			25.40	22.48	4.77	5.10578844	6.30102999	19.41	-3.35032625
A05	HEX	spalt	E5			21.41				6.30102999	17.62	
A06	HEX	spalt	E5			20.63				6.30102999	17.57	intercept
A07	HEX	spalt	E4			25.80	24.69	1.96	4.44521152	5.30102999	25.40	39.5870685
A08	HEX	spalt	E4			24.44				5.30102999	21.41	
A09	HEX	spalt	E4			23.84				5.30102999	20.63	Efficiency
B01	HEX	spalt	66RF	forewing	4	28.01	27.73	OK	3.53816727	4.30102999	25.80	0.99414193
B02	HEX	spalt				27.61				4.30102999	24.44	99%
B03	HEX	spalt				27.58				4.30102999	23.84	
B04	HEX	spalt	63RF	forewing	4	28.10	27.88	OK	3.49436274	3.30102999	28.61	R^2
B05	HEX	spalt				27.77				3.30102999	28.67	0.89577073
B06	HEX	spalt				27.77				3.30102999	28.61	
B07	HEX	spalt	216RF	forewing	5	32.17	30.84	2.50	2.60946775			
B08	HEX	spalt				30.69						
B09	HEX	spalt				29.67						
C01	HEX	spalt	66RH	hindwing	4	27.87	27.31	OK	3.66513619			
C02	HEX	spalt				27.11						
C03	HEX	spalt				26.94						
C04	HEX	spalt	63RH	hindwing	4	28.59	27.56	1.80	3.58869336			
C05	HEX	spalt				27.30						
C06	HEX	spalt				26.80						
C07	HEX	spalt	216RH	hindwing	5	31.28	29.95	2.43	2.87559141			
C08	HEX	spalt				29.73						

Well	Fluor	Target	Sample	tissue	stage	Cq	mean Cq	range	SQ	SQ	Cq	slope
C09	HEX	spalt				28.85						
D01	HEX	spalt	67RF	forewing	4	30.89	30.26	OK	2.78311995			
D02	HEX	spalt				30.10						
D03	HEX	spalt				29.79						
D04	HEX	spalt	217RF	forewing	5	32.18	29.90	4.07	2.89214649			
D05	HEX	spalt				29.41						
D06	HEX	spalt				28.11						
D07	HEX	spalt	212RF	forewing	5	31.28	30.21	2.11	2.7989256			
D08	HEX	spalt				30.19						
D09	HEX	spalt				29.17						
E01	HEX	spalt	67RH	hindwing	4	30.25	30.12	OK	2.82535180			
E02	HEX	spalt				29.99						
E03	HEX	spalt				30.13						
E04	HEX	spalt	217RH	hindwing	5	33.45	30.19	5.25	2.80364544			
E05	HEX	spalt				28.94						
E06	HEX	spalt				28.20						
E07	HEX	spalt	212RH	hindwing	5	31.98	30.10	3.75	2.83067528			
E08	HEX	spalt				30.10						
E09	HEX	spalt				28.23						
F01	HEX	spalt	45RF	forewing	4	28.36	27.77	OK	3.52658928			
F02	HEX	spalt				27.46						
F03	HEX	spalt				27.50						
F04	HEX	spalt	218RF	forewing	5	37.26	33.28	6.06	1.88315859			
F05	HEX	spalt				31.37						
F06	HEX	spalt				31.20						
F07	HEX	spalt	NRT			36.07	33.85	3.83	1.71122271			
F08	HEX	spalt	NRT			33.25						
F09	HEX	spalt	NRT			32.24						
G01	HEX	spalt	45RH	hindwing	4	28.51	28.07	OK	3.43818026			
G02	HEX	spalt				27.82						
G03	HEX	spalt				27.88						

Well	Fluor	Target	Sample	tissue	stage	Cq	mean Cq	range	SQ	SQ	Cq	slope
G04	HEX	spalt	218RH	hindwing	5	37.13	33.84	5.30	1.71578292			
G05	HEX	spalt				32.55						
G06	HEX	spalt				31.83						
G07	HEX	spalt	NTC				#DIV/0!	OK	#DIV/0!			
G08	HEX	spalt	NTC									
G09	HEX	spalt	NTC									
H01	HEX	spalt	E3			28.61	28.63	OK	3.26998949			
H02	HEX	spalt	E3			28.67						
H03	HEX	spalt	E3			28.61						
A01	Texas Red	dpp	E6			19.81	18.59	1.84	6.37964774			
A02	Texas Red	dpp	E6			17.99						
A03	Texas Red	dpp	E6			17.97				SQ	Cq	slope
A04	Texas Red	dpp	E5			25.72	22.80	4.78	5.11271310	6.30102999	19.81	-3.32015375
A05	Texas Red	dpp	E5			21.73				6.30102999	17.99	
A06	Texas Red	dpp	E5			20.94				6.30102999	17.97	intercept
A07	Texas Red	dpp	E4			26.08	25.03	1.88	4.44181054	5.30102999	25.72	39.7733852
A08	Texas Red	dpp	E4			24.80				5.30102999	21.73	
A09	Texas Red	dpp	E4			24.20				5.30102999	20.94	Efficiency
B01	Texas Red	dpp	66RF	forewing	4	24.44	24.23	OK	4.68171628	4.30102999	26.08	1.00037049
B02	Texas Red	dpp				24.10				4.30102999	24.80	100%
B03	Texas Red	dpp				24.14				4.30102999	24.20	
B04	Texas Red	dpp	63RF	forewing	4	22.70	22.53	OK	5.19402079	3.30102999	28.87	R^2
B05	Texas Red	dpp				22.41				3.30102999	29.03	0.89560127
B06	Texas Red	dpp				22.47				3.30102999	28.85	
B07	Texas Red	dpp	216RF	forewing	5	33.05	32.08	1.94	2.31636374			
B08	Texas Red	dpp				32.09						
B09	Texas Red	dpp				31.11						
C01	Texas Red	dpp	66RH	hindwing	4	24.58	23.98	OK	4.75589063			
C02	Texas Red	dpp				23.75						
C03	Texas Red	dpp				23.62						
C04	Texas Red	dpp	63RH	hindwing	4	23.97	22.90	1.83	5.0810706			

Well	Fluor	Target	Sample	tissue	stage	Cq	mean Cq	range	SQ	SQ	Cq	slope
C05	Texas Red	dpp				22.61						
C06	Texas Red	dpp				22.13						
C07	Texas Red	dpp	216RH	hindwing	5	28.50	27.24	2.35	3.77635606			
C08	Texas Red	dpp				27.06						
C09	Texas Red	dpp				26.15						
D01	Texas Red	dpp	67RF	forewing	4	27.27	26.58	OK	3.97378960			
D02	Texas Red	dpp				26.32						
D03	Texas Red	dpp				26.14						
D04	Texas Red	dpp	217RF	forewing	5	27.06	24.67	4.23	4.54953203			
D05	Texas Red	dpp				24.12						
D06	Texas Red	dpp				22.83						
D07	Texas Red	dpp	212RF	forewing	5	26.82	25.67	2.27	4.24710324			
D08	Texas Red	dpp				25.65						
D09	Texas Red	dpp				24.55						
E01	Texas Red	dpp	67RH	hindwing	4	24.67	24.57	OK	4.57764117			
E02	Texas Red	dpp				24.49						
E03	Texas Red	dpp				24.56						
E04	Texas Red	dpp	217RH	hindwing	5	29.61	26.35	5.21	4.04261255			
E05	Texas Red	dpp				25.05						
E06	Texas Red	dpp				24.40						
E07	Texas Red	dpp	212RH	hindwing	5	28.30	26.38	3.83	4.03297436			
E08	Texas Red	dpp				26.39						
E09	Texas Red	dpp				24.46						
F01	Texas Red	dpp	45RF	forewing	4	24.54	23.89	OK	4.78404612			
F02	Texas Red	dpp				23.61						
F03	Texas Red	dpp				23.52						
F04	Texas Red	dpp	218RF	forewing	5	30.61	27.57	4.73	3.67695032			
F05	Texas Red	dpp				26.20						
F06	Texas Red	dpp				25.88						
F07	Texas Red	dpp	NRT			35.90	33.83	3.35	1.78902780			
F08	Texas Red	dpp	NRT			33.05						

Well	Fluor	Target	Sample	tissue	stage	Cq	mean Cq	range	SQ	SQ	Cq	slope
F09	Texas Red	dpp	NRT			32.55						
G01	Texas Red	dpp	45RH	hindwing	4	24.35	23.86	OK	4.79332730			
G02	Texas Red	dpp				23.56						
G03	Texas Red	dpp				23.66						
G04	Texas Red	dpp	218RH	hindwing	5	34.29	31.71	4.33	2.42953923			
G05	Texas Red	dpp				30.87						
G06	Texas Red	dpp				29.96						
G07	Texas Red	dpp	NTC				38.53	OK	0.37390432			
G08	Texas Red	dpp	NTC			38.29						
G09	Texas Red	dpp	NTC			38.78						
H01	Texas Red	dpp	E3			28.87	28.92	OK	3.26994859			
H02	Texas Red	dpp	E3			29.03						
H03	Texas Red	dpp	E3			28.85						

:

v

:

✓

:

✓

;

✓

specimen	tissue	stage	EF1a	optix	dpp	spalt
45	forewing	4	7.08119933	2.61113065	4.78404612	3.52658928
63	forewing	4	7.53959732	2.46389491	5.19402079	3.49436274
66	forewing	4	7.18406609	2.80922007	4.68171628	3.53816727
67	forewing	4	6.38238009	2.09427232	3.97378960	2.78311995
212	forewing	5	6.72314452	2.15559548	4.24710324	2.7989256
217	forewing	5	6.98092905	2.31474567	4.54953203	2.89214649
218	forewing	5	6.31832029	1.56820110	3.67695032	1.88315859
45	hindwing	4	7.15207026	1.77754100	4.79332730	3.43818026
63	hindwing	4	7.40649324	2.98186035	5.0810706	3.58869336
66	hindwing	4	7.12332708	2.39207592	4.75589063	3.66513619
67	hindwing	4	6.96529265	2.06901009	4.57764117	2.82535180
212	hindwing	5	6.46965683	2.67525680	4.03297436	2.83067528
216	hindwing	5	6.46521760	1.90634038	3.77635606	2.87559141
217	hindwing	5	6.38266996	2.11657409	4.04261255	2.80364544
218	hindwing	5	5.08392914	1.41899515	2.42953923	1.71578292
No RT			1.56737823	1.60760334	1.78902780	1.71122271

specimen	tissue	stage	EF1a	optix	dpp	spalt
45	forewing	4	7.08119933	2.61113065	4.78404612	3.52658928
63	forewing	4	7.53959732	2.46389491	5.19402079	3.49436274
66	forewing	4	7.18406609	2.80922007	4.68171628	3.53816727
67	forewing	4	6.38238009	2.09427232	3.97378960	2.78311995
212	forewing	5	6.72314452	2.15559548	4.24710324	2.7989256
217	forewing	5	6.98092905	2.31474567	4.54953203	2.89214649
218	forewing	5	6.31832029	1.56820110	3.67695032	1.88315859
45	hindwing	4	7.15207026	1.77754100	4.79332730	3.43818026
63	hindwing	4	7.40649324	2.98186035	5.0810706	3.58869336
66	hindwing	4	7.12332708	2.39207592	4.75589063	3.66513619
67	hindwing	4	6.96529265	2.06901009	4.57764117	2.82535180
212	hindwing	5	6.46965683	2.67525680	4.03297436	2.83067528
216	hindwing	5	6.46521760	1.90634038	3.77635606	2.87559141
217	hindwing	5	6.38266996	2.11657409	4.04261255	2.80364544
218	hindwing	5	5.08392914	1.41899515	2.42953923	1.71578292